

SECURO DRAIN®

LOW-VACUUM WOUND DRAINAGE

Closed drainage system with manual low-vacuum regulation for a gentle drainage of blood and secretion following e.g. dissection of axillary lymph nodes or mammoplasty

FIELDS OF APPLICATION:

- Abdominal Surgery
- Mamma Surgery
- Plastic Surgery
- Neurosurgery



SECURO DRAIN® silicone reservoir for manual vacuum regulation

FEATURES:

- the non-adhesive silicone material allows an easy, pain-avoiding drainage removal
- the ribbed inner wall design guarantees a problem-free drainage, even when the drainage tube is collapsed
- internal anti-reflux valve prevents backflow of fluid to the patient
- transparent reservoir and graduations make it easy to identify and measure exudate

SECURO DRAIN®

Description	Ref.No.	Size	Length of Perforation	Packaging EA / Box	Drainage-Set with 150cc-Reservoir	Drainage-Set with 400cc-Reservoir	Packaging EA / Box
SECURO DRAIN® Type: flat, 3/4 perforation	01.3307	7 mm	15 cm (3/4)	20	01.3357		5
	01.3310	10 mm	15 cm (3/4)	20	01.3360	01-3360-4	5 / 2
SECURO DRAIN® Typ: flat, complete perforation	01.3404	4 mm	15 cm (4/4)	20	01.3454		5
	01.3407	7 mm	20 cm (4/4)	20	01.3457		5
	01.3410	10 mm	20 cm (4/4)	20	01.3460		5

SECURO DRAIN®, rund

Description	Ref.No.	Size	Length of Perforation	Total Length	Packaging EA / Box
SECURO DRAIN® Type: round	01.3508	7 F	23 cm	125 cm	20
	01.3510	10 F	23 cm	125 cm	20
	01.3515	15 F	23 cm	125 cm	20

SECURO DRAIN® Reservoir

Description	Ref.No.	Contents	Packaging EA / Box
SECURO DRAIN® Reservoir	01.3300	150 cc	5
	01.3303	400 cc	2

STERILE EO LATEX-FREE

SPECIFICATIONS:

- Material:
Drain: Silicone, radiopaque
- Dimensions Flat Drain: 80cm (total)
Flat part: 20cm
Extension line: 60cm